



Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
Cardiac Markers							
CCK-402	CK-MB Rapid Test Cassette	WB / S/ P	Cassette	25 T	5ng/mL	CE	
CCR-401	C-reactive protein Rapid Test Dipstick	WB / S/ P	Dipstick	20 T	10µg/mL	CE	
CCR-402	C-reactive protein Rapid Test Cassette	WB / S/ P	Cassette	20 T	10µg/mL*, 0.5 µg/mL	CE	
CCR-U401	C-reactive protein Semi-Quantitative Rapid Test Dipstick	WB / S/ P	Dipstick	20 T	10µg/mL	CE	
CCR-T402	C-reactive protein Semi-Quantitative Rapid Test Cassette	WB / S/ P	Cassette	20 T	10µg/mL	CE	
CDM-402	D-dimer Rapid Test Cassette	WB / P	Cassette	10 T	500ng/mL	CE	
CFA-402	H-FABP Rapid Test Cassette	WB / S/ P	Cassette	10 T	8ng/mL	CE	
CFC-425	H-FABP and cTnI Combo Rapid Test Cassette	WB / S/ P	Cassette	10 T	See Insert	CE	
CMA-435	Myoglobin/CK-MB/Troponin I Combo Rapid Test Cassette	WB / S/ P	Cassette	25 T	See Insert	CE	
CMY-402	Myoglobin Rapid Test Cassette	WB / S/ P	Cassette	25 T	50ng/mL	CE	
CPC-302	PCT Rapid Test Cassette	S / P	Cassette	10 T	See Insert	CE	
CPC-402	PCT Rapid Test Cassette	WB / S/ P	Cassette	10 T	See Insert	CE	
CTI-402	cTnI Rapid Test Cassette	WB / S/ P	Cassette	25 T	0.5ng/mL	CE	
CBN-402	NT-proBNP Rapid Test Cassette	WB/S/P	Cassette	10 T	0.45ng/mL	CE	
CMA-445	H-FABP and Myoglobin/CK-MB/Cardiac Troponin I Combo Rapid Test Cassette	WB/S/P	Cassette	10 T	See Insert	CE	