



Other

Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
OHS-901	Human Semen Rapid Test Dipstick	Stain	Dipstick	50 T	4ng/mL	CE	
OHS-902	Human Semen Rapid Test Cassette	Stain	Cassette	25 T	4ng/mL	CE	
OMIA-101	Micro-Albumin Semi-Quantitative Rapid Test Dipstick	Urine	Dipstick	50 T	20µg/mL	CE	
OMIA-102	Micro-Albumin Semi-Quantitative Rapid Test Cassette	Urine	Cassette	25T	20µg/mL	CE	
OMAL-101	Micro-Albumin Qualitative Rapid Test Dipstick	Urine	Dipstick	50 T	20µg/mL	CE	
OMAL-102	Micro-Albumin Qualitative Rapid Test Cassette	Urine	Cassette	25 T	20µg/mL	CE	
ORF-302	Rheumatoid Factor Rapid Test Cassette	S/P	Cassette	20T	25IU/mL	CE	
OSC-902	Sperm Concentration Rapid Test Cassette	Sperm	Cassette	2 T	15million/mL	CE	
OSC-902H	Sperm Concentration Rapid Test Cassette	Sperm	Cassette	1 T/2T	15million/mL	Non-CE	
OSP-902	SP-10 Male Fertility Rapid Test Cassette	Sperm	Cassette	1 T/2T	15 million/mL	CE	
OSA-402	SAA Rapid Test Cassette	WB/S/P	Cassette	10 T	10µg/mL	CE	
OCSC-425	SAA and CRP Combo Rapid Test Cassette	WB/S/P	Cassette	10 T	10µg/mL	CE	
OTS-402	TSH Rapid Test Cassette	WB / S/ P	Cassette	25 T	5µIU/mL	CE	
OTS-402H	TSH Rapid Test Cassette	WB	Cassette	1 T	5µIU/mL	Non-CE	
OVD-402	Vitamin D Rapid Test Cassette	WB	Cassette	10 T	30ng/mL	CE	
OHA-402	HbA1c Rapid Test Cassette	WB	Cassette	10 T	6%~9%	CE	
U031-04H	Urinary Tract Infections Test	Urine	Dipstick	1T	Analytes: LEU, BLO,NIT,PRO) Analytes: (ASC, BLO, BIL, URO,KET,GLU, PRO, NIT, LEU, pH, SG,ALB,CRE)	Non-CE	
U031-01	Urinalysis Reagent Dipsticks - 1 parameter	Urine	Dipstick	100 T		CE	
U031-02	Urinalysis Reagent Dipsticks - 2 parameters	Urine	Dipstick	100 T		CE	
U031-03	Urinalysis Reagent Dipsticks - 3 parameters	Urine	Dipstick	100 T		CE	
U031-04	Urinalysis Reagent Dipsticks - 4 parameters	Urine	Dipstick	100 T		CE	
U031-05	Urinalysis Reagent Dipsticks - 5 parameters	Urine	Dipstick	100 T		CE	
U031-06	Urinalysis Reagent Dipsticks - 6 parameters	Urine	Dipstick	100 T		CE	
U031-07	Urinalysis Reagent Dipsticks - 7 parameters	Urine	Dipstick	100 T		CE	
U031-08	Urinalysis Reagent Dipsticks - 8 parameters	Urine	Dipstick	100 T		CE	
U031-09	Urinalysis Reagent Dipsticks - 9 parameters	Urine	Dipstick	100 T		CE	
U031-10	Urinalysis Reagent Dipsticks - 10 parameters	Urine	Dipstick	100 T		CE	
U031-11	Urinalysis Reagent Dipsticks - 11 parameters	Urine	Dipstick	100 T		CE	

Other

OAL-914	Breast milk Alcohol Rapid Test Panel	Breast Milk	Panel	25 T	0.01%	CE	
OBS-901	Blood Stain Rapid Test Dipstick	Stain	Dipstick	50T	50ng/mL	CE	
OBS-902	Blood Stain Rapid Test Cassette	Stain	Cassette	25 T	50ng/mL	CE	
OCAL-602	Calprotectin test cassette	Feces	Cassette	10 T	50µg/g	CE	
OLF-602	Lactoferrin Rapid Test Cassette	Feces	Cassette	10 T	100 ng/mL	CE	
OIGE-402	IgE Rapid Test Cassette	WB/S/P	Cassette	10 T	200 IU/mL	CE	
OFE-T402	Ferritin Semi-Quantitative Rapid Test	WB / S/ P	Cassette	10 T	See Insert	CE	
OFE-402	Ferritin Rapid test cassette	WB / S/ P	Cassette	10 T	30ng/mL	CE	
OFE-402H	Ferritin Rapid test cassette	WB	Cassette	1 T	30ng/mL	Non-CE	