





Tumor Markers

Tumor Markers							
TAF-402	AFP Rapid Test Cassette	WB / S/ P	Cassette	25 T	10ng/mL	CE	
TCE-402	CEA Rapid Test Cassette	WB / S/ P	Cassette	25 T	5ng/mL	CE	
TFO-602	FOB Rapid Test Cassette	Feces	Cassette	25 T	40ng/mL,50ng/mL*	CE	
TFO-602H	FOB Rapid Test Cassette	Feces	Cassette	1 T		Non-CE	
TFO-601	FOB Rapid Test Dipstick	Feces	Dipstick	25 T		CE	
TFO-601H	FOB Rapid Test Dipstick	Feces	Dipstick	1 T		Non-CE	
TPS-402	PSA Rapid Test Cassette	WB / S/ P	Cassette	25 T	3ng/mL	Non-CE	
TPS-Q402	PSA Qualitative Rapid Test Cassette	WB / S/ P	Cassette	25 T	4ng/mL	Non-CE	
TPS-Q401	PSA Qualitative Rapid Test Dipstick	WB / S/ P	Dipstick	25 T	4ng/mL	Non-CE	
T125-402	CA125 Rapid Test Cassette	WB/S/P	Cassette	10 T	40U/mL	CE	
T153-402	CA15-3 Rapid Test Cassette	WB/S/P	Cassette	10 T	30U/mL	CE	
T199-402	CA19-9 Rapid Test Cassette	WB/S/P	Cassette	10 T	40U/mL	CE	
TTFC-625	Transferrin and FOB Combo Rapid Test Cassette	Feces	Cassette	25 T	See Insert	CE	
THPB-625	Hb+Hb-Hp Combo Rapid Test Cassette	Feces	Cassette	25 T	50ng/mL	CE	
TTFH-635	Transferrin/FOB and Hb-Hp Combo Rapid Test Cassette	Feces	Cassette	25 T	See Insert	CE	