
























Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture		
<b>Drugs of Abuse</b>									
DABP-102	AB-PINACA (ABP) Rapid Test Cassette	Urine	Cassette	25 T	10 ng/mL*	CE	 <p>Cassette</p>		
DABP-114	AB-PINACA (ABP) Rapid Test Panel	Urine	Panel	25 T		CE			
DABP-101	AB-PINACA (ABP) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DAC-102	Acetaminophen (ACE) Rapid Test Cassette	Urine	Cassette	25 T	5000 ng/mL*	CE		 <p>Cassette</p>	
DAC-114	Acetaminophen (ACE) Rapid Test Panel	Urine	Panel	25 T		CE			
DAC-101	Acetaminophen (ACE) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DACL-102	7-Aminoclonazepam (7-ACL) Rapid Test Cassette	Urine	Cassette	25 T	300 ng/mL*	CE			 <p>Cassette</p>
DACL-114	7-Aminoclonazepam (7-ACL) Rapid Test Panel	Urine	Panel	25 T		CE			
DACL-101	7-Aminoclonazepam (7-ACL) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DAL-102	Alcohol (ALC) Rapid Test Cassette	Urine	Cassette	25 T	0.02%	CE	 <p>Cassette</p>		
DAL-114	Alcohol (ALC) Rapid Test Panel	Urine	Panel	25 T		CE			
DAL-101	Alcohol (ALC) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DAM-102	Amphetamine (AMP) Rapid Test Cassette	Urine	Cassette	25 T	300 ng/mL, 500 ng/mL, 1000 ng/mL*	CE/FDA		 <p>Cassette</p>	
DAM-114	Amphetamine (AMP) Rapid Test Panel	Urine	Panel	25 T		CE/FDA			
DAM-101	Amphetamine (AMP) Rapid Test Dipstick	Urine	Dipstick	50 T		CE/FDA			
DAP-102	α-Pyrrolidinovalerophenone (α -PVP ) Rapid	Urine	Cassette	25 T	1000 ng/mL*, 500ng/mL, 300ng/mL	CE			 <p>Cassette</p>
DAP-114	α-Pyrrolidinovalerophenone (α -PVP ) Rapid T	Urine	Panel	25 T		CE			
DAP-101	α-Pyrrolidinovalerophenone (α -PVP ) Rapid T	Urine	Dipstick	50 T		CE			
DBA-102	Barbiturate (BAR) Rapid Test Cassette	Urine	Cassette	25 T	200 ng/mL, 300 ng/mL*	CE	 <p>Cassette</p>		
DBA-114	Barbiturate (BAR) Rapid Test Panel	Urine	Panel	25 T		CE			
DBA-101	Barbiturate (BAR) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DBU-102	Buprenorphine (BUP) Rapid Test Cassette	Urine	Cassette	25 T	10 ng/mL*	CE		 <p>Cassette</p>	
DBU-114	Buprenorphine (BUP) Rapid Test Panel	Urine	Panel	25 T		CE			
DBU-101	Buprenorphine (BUP) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DBZ-102	Benzodiazepines (BZO) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL, 200 ng/mL, 300 ng/mL*	CE/FDA			 <p>Cassette</p>
DBZ-114	Benzodiazepines (BZO) Rapid Test Panel	Urine	Panel	25 T		CE/FDA			
DBZ-101	Benzodiazepines (BZO) Rapid Test Dipstick	Urine	Dipstick	50 T		CE/FDA			
DCA-102	Cathine (CAT) Rapid Test Cassette	Urine	Cassette	25 T	150 ng/mL*	CE	 <p>Cassette</p>		
DCA-114	Cathine (CAT) Rapid Test Panel	Urine	Panel	25 T		CE			
DCA-101	Cathine (CAT) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DCNB-102	Cannabinol (CNB) Rapid Test Cassette	Urine	Cassette	25 T	500 ng/mL	CE		 <p>Cassette</p>	
DCNB-101	Cannabinol (CNB) Rapid Test Panel	Urine	Panel	25 T		CE			
DCNB-114	Cannabinol (CNB) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DCAF-102	Caffeine(CAF) Rapid Test Cassette	Urine	Cassette	25 T	1000 ng/mL*	CE			 <p>Cassette</p>
DCAF-114	Caffeine(CAF) Rapid Test Panel	Urine	Panel	25 T		CE			
DCAF-101	Caffeine(CAF) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			



Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture		
DCL-102	Clonazepam (CLO) Rapid Test Cassette	Urine	Cassette	25 T	400 ng/mL*	CE	 Panel		
DCL-114	Clonazepam (CLO) Rapid Test Panel	Urine	Panel	25 T		CE			
DCL-101	Clonazepam (CLO) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DCO-102	COCAINE (COC) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL, 150 ng/mL, 300 ng/mL*	CE/FDA			
DCO-114	COCAINE (COC) Rapid Test Panel	Urine	Panel	25 T		CE/FDA			
DCO-101	COCAINE (COC) Rapid Test Dipstick	Urine	Dipstick	50 T		CE/FDA			
DCT-102	Cotinine (COT) Rapid Test Cassette	Urine	Cassette	25 T	10 ng/mL,50 ng/mL, 100 ng/mL, 200 ng/mL*, 500 ng/mL	CE			
DCT-114	Cotinine (COT) Rapid Test Panel	Urine	Panel	25 T		CE			
DCT-101	Cotinine (COT) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DDI-102	Diazepam (DIA) Rapid Test Cassette	Urine	Cassette	25 T	300 ng/mL*	CE			
DDI-114	Diazepam (DIA) Rapid Test Panel	Urine	Panel	25 T		CE			
DDI-101	Diazepam (DIA) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DED-102	Ethylenediamine-dimethylphosphinic acid (EDDP) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL*	CE			
DED-114	Ethylenediamine-dimethylphosphinic acid (EDDP) Rapid Test Panel	Urine	Panel	25 T		CE			
DED-101	Ethylenediamine-dimethylphosphinic acid (EDDP) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DET-102	Ethyl Glucuronide (ETG) Rapid Test Cassette	Urine	Cassette	25 T	500 ng/mL*, 1000 ng/mL	CE			 Cassette
DET-114	Ethyl Glucuronide (ETG) Rapid Test Panel	Urine	Panel	25 T		CE			
DET-101	Ethyl Glucuronide (ETG) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DFY-102	Fentanyl (FYL) Rapid Test Cassette	Urine	Cassette	25 T	10 ng/mL ,20 ng/mL*	CE			
DFY-114	Fentanyl (FYL) Rapid Test Panel	Urine	Panel	25 T		CE			
DFY-101	Fentanyl (FYL) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DGAB-102	Gabapentin(GAB) Rapid Test Cassette	Urine	Cassette	25 T	2000 ng/mL	CE			
DGAB-114	Gabapentin(GAB) Rapid Test Panel	Urine	Panel	25 T		CE			
DGAB-101	Gabapentin(GAB) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DGHB-102	γ-Hydroxybutyric acid (GHB) Rapid Test Cassette	Urine	Cassette	25 T	10µg/mL	CE			
DKE-102	Ketamine (KET)Rapid Test Cassette	Urine	Cassette	25 T	500 ng/mL, 1000 ng/mL*	CE			
DKE-114	Ketamine (KET)Rapid Test Panel	Urine	Panel	25 T		CE			
DKE-101	Ketamine (KET)Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DMA-102	6-Monoacetylmorphine (6-MAM) Rapid Test Cassette	Urine	Cassette	25 T	10 ng/mL*	CE			
DMA-114	6-Monoacetylmorphine (6-MAM)Rapid Test Panel	Urine	Panel	25 T		CE			
DMA-101	6-Monoacetylmorphine(6-MAM) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DMCA-102	Methcathinone(MCAT) Rapid Test Cassette	Urine	Cassette	25 T	500 ng/mL*	CE		 Dipstick	
DMCA-114	Methcathinone(MCAT) Rapid Test Panel	Urine	Panel	25 T		CE			
DMCA-101	Methcathinone(MCAT) Rapid Test Dipstick	Urine	Dipstick	50 T		CE			
DMD-102	Ecstasy (MDMA) Rapid Test Cassette	Urine	Cassette	25 T	500 ng/mL*, 1000 ng/mL	CE/FDA			
DMD-114	Ecstasy (MDMA) Rapid Test Panel	Urine	Panel	25 T		CE/FDA			
DMD-101	Ecstasy (MDMA) Rapid Test Dipstick	Urine	Dipstick	50 T		CE/FDA			

Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
DMDA-102	Tenamfetamine (MDA) Rapid Test Cassette	Urine	Cassette	25 T	500 ng/mL*	CE	 Panel
DMDA-114	Tenamfetamine (MDA) Rapid Test Panel	Urine	Panel	25 T		CE	
DMDA-101	Tenamfetamine (MDA) Rapid Test Dipstick	Urine	Dipstick	50 T		CE	
DMDP-102	Methylenedioxypropyvalerone (MDPV) Rapid Test Cassette	Urine	Cassette	25 T	1000 ng/mL*, 500ng/mL, 300ng/mL	CE	
DMDP-114	Methylenedioxypropyvalerone (MDPV) Rapid Test Panel	Urine	Panel	25 T		CE	
DMDP-101	Methylenedioxypropyvalerone (MDPV) Rapid Test Dipstick	Urine	Dipstick	50 T		CE	
DMP-102	Methylphenidate(MPD) Rapid Test Cassette	Urine	Cassette	25 T	1000 ng/mL* 300 ng/mL	CE	
DMP-114	Methylphenidate(MPD) Rapid Test Panel	Urine	Panel	25 T		CE	
DMP-101	Methylphenidate(MPD) Rapid Test Dipstick	Urine	Dipstick	50 T		CE	
DMPR-102	Meperidine (MPRD) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL	CE	
DMPR-114	Meperidine (MPRD) Rapid Test Panel	Urine	Panel	25 T		CE	
DMPR-101	Meperidine (MPRD) Rapid Test Dipstick	Urine	Dipstick	50 T		CE	
DME-102	Methamphetamine (MET) Rapid Test Cassette	Urine	Cassette	25 T	300 ng/mL, 500 ng/mL, 1000 ng/mL*	CE/FDA	
DME-114	Methamphetamine (MET) Rapid Test Panel	Urine	Panel	25 T		CE/FDA	
DME-101	Methamphetamine (MET) Rapid Test Dipstick	Urine	Dipstick	50 T		CE/FDA	
DMEP-102	Mephedrone HCl(MEP) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL	CE	
DMEP-114	Mephedrone HCl(MEP) Rapid Test Panel	Urine	Panel	25 T		CE	
DMEP-101	Mephedrone HCl(MEP) Rapid Test Dipstick	Urine	Dipstick	50 T		CE	
DMES-102	Mescaline (MES) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL*, 300ng/mL	CE	
DMES-114	Mescaline (MES) Rapid Test Panel	Urine	Panel	25 T		CE	
DMES-114	Mescaline (MES) Rapid Test Dipstick	Urine	Dipstick	50 T		CE	
<b>Drugs of Abuse</b>							
DMO-102	Morphine (MOP) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL, 300 ng/mL*	CE	 Cassette
DMO-114	Morphine (MOP) Rapid Test Panel	Urine	Panel	25 T		CE	
DMO-101	Morphine (MOP) Rapid Test Dipstick	Urine	Dipstick	50 T		CE	
DMQ-102	Methaqualone (MQL) Rapid Test Cassette	Urine	Cassette	25 T	300 ng/mL*	CE	
DMQ-114	Methaqualone (MQL) Rapid Test Panel	Urine	Panel	25 T		CE	
DMQ-101	Methaqualone (MQL) Rapid Test Dipstick	Urine	Dipstick	50 T		CE	
DMT-102	Methadone (MTD) Rapid Test Cassette	Urine	Cassette	25 T	200 ng/mL, 300 ng/mL*	CE	
DMT-114	Methadone (MTD) Rapid Test Panel	Urine	Panel	25 T		CE	
DMT-101	Methadone (MTD) Rapid Test Dipstick	Urine	Dipstick	50 T		CE	
DOP-102	Opiates (OPI) Rapid Test Cassette	Urine	Cassette	25 T	2000 ng/mL* 300 ng/mL	CE/FDA	
DOP-114	Opiates (OPI) Rapid Test Panel	Urine	Panel	25 T		CE/FDA	
DOP-101	Opiates (OPI) Rapid Test Dipstick	Urine	Dipstick	50 T		CE/FDA	

Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture	
DOX-102	Oxycodone (OXY) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL*	CE	 <p>Dipstick</p>	
DOX-114	Oxycodone (OXY) Rapid Test Panel	Urine	Panel	25 T		CE		
DOX-101	Oxycodone (OXY) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DPC-102	Phencyclidine (PCP) Rapid Test Cassette	Urine	Cassette	25 T	25 ng/mL*	CE/FDA		
DPC-114	Phencyclidine (PCP) PCP Rapid Test Panel	Urine	Panel	25 T		CE/FDA		
DPC-101	Phencyclidine (PCP) PCP Rapid Test Dipstick	Urine	Dipstick	50 T		CE/FDA		
DPG-102	Pregabalin (PGB) Rapid Test Cassette	Urine	Cassette	25 T	50µg/mL*	CE		
DPG-114	Pregabalin (PGB) Rapid Test Panel	Urine	Panel	25 T		CE		
DPG-101	Pregabalin (PGB) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DPP-102	Propoxyphene (PPX) Rapid Test Cassette	Urine	Cassette	25 T	300 ng/mL*	CE		
DPP-114	Propoxyphene (PPX) Rapid Test Panel	Urine	Panel	25 T		CE		
DPP-101	Propoxyphene (PPX) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DSM-102	Synthetic Marijuana (K2) Rapid Test Cassette	Urine	Cassette	25 T	50 ng/mL*, 30 ng/mL	CE		
DSM-114	Synthetic Marijuana (K2) Rapid Test Panel	Urine	Panel	25 T		CE		
DSM-101	Synthetic Marijuana (K2) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DTC-102	Tricyclic Antidepressants (TCA) Rapid Test Cassette	Urine	Cassette	25 T	500 ng/mL,1000 ng/mL*	CE	 <p>Panel</p>	
DTC-114	Tricyclic Antidepressants (TCA) Rapid Test Panel	Urine	Panel	25 T		CE		
DTC-101	Tricyclic Antidepressants (TCA) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DTH-102	Marijuana (THC) Rapid Test Cassette	Urine	Cassette	25 T	25 ng/mL,20 ng/mL, 50 ng/mL*, 150 ng/mL, 300 ng/mL	CE/FDA		
DTH-114	Marijuana (THC) Rapid Test Panel	Urine	Panel	25 T		CE/FDA		
DTH-101	Marijuana (THC) Rapid Test Dipstick	Urine	Dipstick	50 T		CE/FDA		
DTM-102	Tramadol (TML) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL*,200 ng/mL, 300 ng/mL	CE		
DTM-114	Tramadol (TML) Rapid Test Panel	Urine	Panel	25 T		CE		
DTM-101	Tramadol (TML) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DTRO-102	Tropicamide(TRO)Rapid Test Cassette	Urine	Cassette	25 T	350 ng/mL*	CE		
DTRO-114	Tropicamide(TRO)Rapid Test Panel	Urine	Panel	25 T		CE		
DTRO-101	Tropicamide(TRO) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DUR-102	UR-144 Rapid Test Cassette	Urine	Cassette	40 T	25 ng/mL	CE		
DUR-114	UR-144 Rapid Test Panel	Urine	Panel	40 T		CE		
DUR-101	UR-144 Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DTZD-102	Trazodone(TZD) Test Cassette	Urine	Cassette	25 T	200 ng/mL	CE		
DTZD-114	Trazodone(TZD) apid Test Panel	Urine	Panel	25 T		CE		
DTZD-101	Trazodone(TZD) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DZAL-102	Zaleplon(ZAL) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL	CE		
DZAL-114	Zaleplon(ZAL) Rapid Test Panel	Urine	Panel	25 T		CE		
DZAL-101	Zaleplon(ZAL) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		

Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture	
DLS-102	Lysergic Acid Diethylamide (LSD) Rapid Test Cassette	Urine	Cassette	25 T	10 ng/mL, 20 ng/mL*, 50 ng/mL	CE	 <p>Dipstick</p>	
DLS-114	Lysergic Acid Diethylamide (LSD) Rapid Test Panel	Urine	Panel	25 T		CE		
DLS-101	Lysergic Acid Diethylamide (LSD) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DZO-102	Zolpidem(ZOL) Rapid Test Cassette	Urine	Cassette	25 T	50 ng/mL*	CE		
DZO-114	Zolpidem(ZOL) Rapid Test Panel	Urine	Panel	25 T		CE		
DZO-101	Zolpidem(ZOL) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DALP-102	Alprazolam(ALP) Rapid Test Cassette	Urine	Cassette	25 T	100 ng/mL*	CE		 <p>Panel</p>
DALP-114	Alprazolam(ALP) Rapid Test Panel	Urine	Panel	25 T		CE		
DALP-101	Alprazolam(ALP) Rapid Test Dipstick	Urine	Dipstick	50 T		CE		
DZOP-101	Zopiclone(ZOP) Rapid Test Dipstick	Urine	Dipstick	50 T	50 ng/mL*	CE		
DZOP-102	Zopiclone(ZOP) Rapid Test Cassette	Urine	Cassette	25 T		CE		
DZOP-114	Zopiclone(ZOP) Rapid Test Panel	Urine	Panel	25 T		CE		
DCFY-101	Carfentanyl(CFYL) Rapid Test Dipstick	Urine	Dipstick	50 T	500 ng/mL*	CE		
DCFY-102	Carfentanyl(CFYL) Rapid Test Cassette	Urine	Cassette	25 T		CE		
DCFY-114	Carfentanyl(CFYL) Rapid Test Panel	Urine	Panel	25 T		CE		
DAT-111	Adulteration Dipstick	Urine	Dipstick	25 T	/	CE		
DOA-125	Multi-Drug 2 Drugs Rapid Test Cassette	Urine	Cassette	25 T	/	CE		
DOA-135	Multi-Drug 3 Drugs Rapid Test Cassette	Urine	Cassette	25 T				
DOA-135	Multi-Drug 3 Drugs Rapid Test Cassette(α-PVP, MDPV, MCAT)	Urine	Cassette	25 T				
DOA-135	Multi-Drug 3 Drugs Rapid Test Cassette(ABP, UR-144, JWH)	Urine	Cassette	25 T				
DOA-145	Multi-Drug 4 Drugs Rapid Test Cassette	Urine	Cassette	25 T				
DOA-145	Multi-Drug 4 Drugs Rapid Test Cassette(K2, α-PVP, MDPV, MCAT)	Urine	Cassette	25 T				
DOA-145	Multi-Drug 4 Drugs Rapid Test Cassette(K2, ABP, UR-144, JWH)	Urine	Cassette	25 T				
DOA-155	Multi-Drug 5 Drugs Rapid Test Cassette	Urine	Cassette	25 T				
DOA-165	Multi-Drug 6 Drugs Rapid Test Cassette	Urine	Cassette	25 T				
DOA-175	Multi-Drug 7 Drugs Rapid Test Cassette	Urine	Cassette	25 T				
DOA-185	Multi-Drug 8 Drugs Rapid Test Cassette	Urine	Cassette	25 T				
DOA-195	Multi-Drug 9 Drugs Rapid Test Cassette	Urine	Cassette	25 T				
DOA-1105	Multi-Drug 10 Drugs Rapid Test Cassette	Urine	Cassette	25 T				
DOA-1115	Multi-Drug 11 Drugs Rapid Test Cassette	Urine	Cassette	25 T				
DOA-1125	Multi-Drug 12 Drugs Rapid Test Cassette	Urine	Cassette	25 T				



Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
DOA-124	Multi-Drug 2 Drugs Rapid Test Panel	Urine	Panel	25 T	/	CE	
DOA-134	Multi-Drug 3 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-134	Multi-Drug 3 Drugs Rapid Test Panel( $\alpha$ -PVP, MDPV, MCAT)	Urine	Panel	25 T			
DOA-134	Multi-Drug 3 Drugs Rapid Test Panel(ABP, UR-144, JWH)	Urine	Panel	25 T			
DOA-144	Multi-Drug 4 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-144	Multi-Drug 4 Drugs Rapid Test Panel(K2, $\alpha$ -PVP, MDPV, MCAT)	Urine	Panel	25 T			
DOA-144	Multi-Drug 4 Drugs Rapid Test Panel(K2,ABP, UR-144, JWH)	Urine	Panel	25 T			
DOA-154	Multi-Drug 5 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-164	Multi-Drug 6 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-174	Multi-Drug 7 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-184	Multi-Drug 8 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-194	Multi-Drug 9 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-1104	Multi-Drug 10 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-1114	Multi-Drug 11 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-1124	Multi-Drug 12 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-1134	Multi-Drug 13 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-1144	Multi-Drug 14 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-1154	Multi-Drug 15 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-1164	Multi-Drug 16 Drugs Rapid Test Panel	Urine	Panel	25 T			
DOA-127-A1	Multi-Drug 2 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T	/	CE/FDA	
DOA-137-A1	Multi-Drug 3 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-147-A1	Multi-Drug 4 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-157-A1	Multi-Drug 5 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-167-A1	Multi-Drug 6 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-177-A1	Multi-Drug 7 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-187-A1	Multi-Drug 8 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-197-A1	Multi-Drug 9 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1107-A1	Multi-Drug 10 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1117-A1	Multi-Drug 11 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1127-A1	Multi-Drug 12 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1137-A1	Multi-Drug 13 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1147-A1	Multi-Drug 14 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1157-A1	Multi-Drug 15 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			


Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
DOA-127-A2	Multi-Drug 2 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T	/	CE	
DOA-137-A2	Multi-Drug 3 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-147-A2	Multi-Drug 4 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-157-A2	Multi-Drug 5 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-167-A2	Multi-Drug 6 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-177-A2	Multi-Drug 7 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-187-A2	Multi-Drug 8 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-197-A2	Multi-Drug 9 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1107-A2	Multi-Drug 10 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1117-A2	Multi-Drug 11 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1127-A2	Multi-Drug 12 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-127-A3	Multi-Drug 2 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-137-A3	Multi-Drug 3 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-147-A3	Multi-Drug 4 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-157-A3	Multi-Drug 5 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-167-A3	Multi-Drug 6 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-177-A3	Multi-Drug 7 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-187-A3	Multi-Drug 8 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-197-A3	Multi-Drug 9 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1107-A3	Multi-Drug 10 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1117-A3	Multi-Drug 11 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1127-A3	Multi-Drug 12 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1137-A3	Multi-Drug 13 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1147-A3	Multi-Drug 14 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-1157-A3	Multi-Drug 15 Drugs Rapid Test 1-Step Cup	Urine	Cup	25 T			
DOA-127-B1	Multi-Drug 2 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T	/	CE	
DOA-137-B1	Multi-Drug 3 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-147-B1	Multi-Drug 4 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-157-B1	Multi-Drug 5 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-167-B1	Multi-Drug 6 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-177-B1	Multi-Drug 7 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-187-B1	Multi-Drug 8 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-197-B1	Multi-Drug 9 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1107-B1	Multi-Drug 10 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1117-B1	Multi-Drug 11 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1127-B1	Multi-Drug 12 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			

Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
DOA-127-B2	Multi-Drug 2 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T	/	CE	
DOA-137-B2	Multi-Drug 3 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-147-B2	Multi-Drug 4 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-157-B2	Multi-Drug 5 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-167-B2	Multi-Drug 6 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-177-B2	Multi-Drug 7 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-187-B2	Multi-Drug 8 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-197-B2	Multi-Drug 9 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1107-B2	Multi-Drug 10 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1117-B2	Multi-Drug 11 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1127-B2	Multi-Drug 12 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1137-B2	Multi-Drug 13 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1147-B2	Multi-Drug 14 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1157-B2	Multi-Drug 15 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DOA-1167-B2	Multi-Drug 16 Drugs Rapid Test 2-Step Cup	Urine	Cup	25 T			
DAL-801	Alcohol Rapid Test Dipstick	Oral Fluid	Dipstick	50 T			
DAL-802	Alcohol Rapid Test Cassette	Oral Fluid	Cassette	25 T	0.02%	CE	
DAL-902	Breath Alcohol Test (With/Without blow bag)	Breath	Cassette	15/20 T	0.08%	CE	
DAP-802	$\alpha$ -Pyrrolidinovalerophenone ( $\alpha$ -PVP) Rapid Test Cassette	Oral Fluid	Cassette	25 T	300 ng/mL	CE	
DAP-803	$\alpha$ -Pyrrolidinovalerophenone ( $\alpha$ -PVP) Rapid Test Midstream	Oral Fluid	Midstream	25 T	300 ng/mL	CE	
DAM-802	Amphetamine (AMP) Rapid Test Cassette	Oral Fluid	Cassette	25 T	50 ng/mL	CE	
DAM-803	Amphetamine (AMP) Rapid Test Midstream	Oral Fluid	Midstream	25 T	50 ng/mL	CE	
DBA-802	Barbiturates (BAR) Rapid Test Cassette	Oral Fluid	Cassette	25 T	50 ng/mL*	CE	
DBA-803	Barbiturates (BAR) Rapid Test Midstream	Oral Fluid	Midstream	25 T	50 ng/mL*	CE	
DBU-802	Buprenorphine (BUP) Rapid Test Cassette	Oral Fluid	Cassette	25 T	10 ng/mL*	CE	
DBU-803	Buprenorphine (BUP) Rapid Test Midstream	Oral Fluid	Midstream	25 T	10 ng/mL*	CE	
DBZ-802	Benzodiazepines (BZO) Rapid Test Cassette	Oral Fluid	Cassette	25 T	20 ng/mL, 30 ng/mL, 50 ng/mL*	CE	
DBZ-803	Benzodiazepines (BZO) Rapid Test Midstream	Oral Fluid	Midstream	25 T	30 ng/mL, 50 ng/mL*	CE	
DCFY-802	Carfentanyl(CFYL) Rapid Test Cassette	Oral Fluid	Cassette	25 T	50 ng/mL	CE	
DCFY-803	Carfentanyl(CFYL) Rapid Test Midstream	Oral Fluid	Midstream	25 T	50 ng/mL	CE	
DCO-802	COCAINE (COC) Rapid Test Cassette	Oral Fluid	Cassette	25 T	20 ng/mL	CE	 <p>Cassette</p>
DCO-803	COCAINE (COC) Rapid Test Midstream	Oral Fluid	Midstream	25 T	20 ng/mL	CE	
DCT-802	Cotinine (COT) Rapid Test Cassette	Oral Fluid	Cassette	25 T	20 ng/mL	CE	
DCT-803	Cotinine (COT) Rapid Test Midstream	Oral Fluid	Midstream	25 T	20 ng/mL	CE	
DFY-802	Fentanyl (FYL) Rapid Test Cassette	Oral Fluid	Cassette	25 T	50ng/ml*, 20ng/ml	CE	
DFY-803	Fentanyl (FYL) Rapid Test Midstream	Oral Fluid	Midstream	25 T	50ng/ml*, 20ng/ml	CE	



Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
DKE-802	Ketamine(KET) Rapid Test Cassette	Oral Fluid	Cassette	25 T	50 ng/mL*	CE	 <p>Midstream</p>
DKE-803	Ketamine(KET) Rapid Test Midstream	Oral Fluid	Midstream	25 T	50 ng/mL*	CE	
DMA-802	6-Monoacetylmorphine (6-MAM) Rapid Test Cassette	Oral Fluid	Cassette	25 T	10 ng/mL	CE	
DMA-803	6-Monoacetylmorphine (7-MAM) Rapid Test Midstream	Oral Fluid	Midstream	25 T	10 ng/mL	CE	
DMD-802	Ecstasy (MDMA) Rapid Test Cassette	Oral Fluid	Cassette	25 T	50 ng/mL	CE	
DMD-803	Ecstasy (MDMA) Rapid Test Midstream	Oral Fluid	Midstream	25 T	50 ng/mL	CE	
DMDP-802	Methylenedioxypropylvalerone (MDPV) Rapid Test Cassette	Oral Fluid	Cassette	25 T	300 ng/mL	CE	
DMDP-803	Methylenedioxypropylvalerone (MDPV) Rapid Test Midstream	Oral Fluid	Midstream	25 T	300 ng/mL	CE	
DME-802	Methamphetamine (MET) Rapid Test Cassette	Oral Fluid	Cassette	25 T	50 ng/mL	CE	
DME-803	Methamphetamine (MET) Rapid Test Midstream	Oral Fluid	Midstream	25 T	50 ng/mL	CE	
DMT-802	Methadone (MTD) Rapid Test Cassette	Oral Fluid	Cassette	25 T	30 ng/mL	CE	
DMT-803	Methadone (MTD) Rapid Test Midstream	Oral Fluid	Midstream	25 T	30 ng/mL	CE	
DOP-802	Opiates (OPI) Rapid Test Cassette	Oral Fluid	Cassette	25 T	40 ng/mL	CE	
DOP-803	Opiates (OPI) Rapid Test Midstream	Oral Fluid	Midstream	25 T	40 ng/mL	CE	
DOX-802	Oxycodone (OXY) Rapid Test Cassette	Oral Fluid	Cassette	25 T	20 ng/mL	CE	
DOX-803	Oxycodone (OXY) Rapid Test Midstream	Oral Fluid	Midstream	25 T	20 ng/mL	CE	
DPC-802	Phencyclidine (PCP) Rapid Test Cassette	Oral Fluid	Cassette	25 T	10 ng/mL	CE	
DPC-803	Phencyclidine (PCP) Rapid Test Midstream	Oral Fluid	Midstream	25 T	10 ng/mL	CE	
DSM-802	Synthetic Marijuana (K2) Rapid Test Cassette	Oral Fluid	Cassette	25 T	25 ng/mL	CE	
DSM-803	Synthetic Marijuana (K2) Rapid Test Midstream	Oral Fluid	Midstream	25 T	25 ng/mL	CE	
DTH-802	Marijuana (THC) Rapid Test Cassette	Oral Fluid	Cassette	25 T	12 ng/mL, 20 ng/mL*	CE	
DTH-803	Marijuana (THC) Rapid Test Midstream	Oral Fluid	Midstream	25 T	50 ng/mL	CE	
DTM-802	Tramadol(TML) Rapid Test Cassette	Oral Fluid	Cassette	25 T	30 ng/mL*	CE	
DTM-803	Tramadol(TML) Rapid Test Midstream	Oral Fluid	Midstream	25 T	30 ng/mL*	CE	
DSD-823-A	Multi-Drug 2 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-833-A	Multi-Drug 3 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-843-A	Multi-Drug 4 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-853-A	Multi-Drug 5 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-863-A	Multi-Drug 6 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-873-A	Multi-Drug 7 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-823-B	Multi-Drug 2 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T	/	CE	
DSD-833-B	Multi-Drug 3 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-843-B	Multi-Drug 4 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-853-B	Multi-Drug 5 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-863-B	Multi-Drug 6 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-873-B	Multi-Drug 7 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	
DSD-883-B	Multi-Drug 8 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T		CE	

Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
DSD-893-B	Multi-Drug 9 Drugs Rapid Test Midstream	Oral Fluid	Midstream	25 T	/	CE	
DSD-825-A	Oral Fluid Multi-Drug 2 Drugs Test A Cassette	Oral Fluid	Cassette	25 T	/	CE	
DSD-835-A	Oral Fluid Multi-Drug 3 Drugs Test A Cassette	Oral Fluid	Cassette	25 T	/	CE	
DSD-845-A	Oral Fluid Multi-Drug 4 Drugs Test A Cassette	Oral Fluid	Cassette	25 T	/	CE	
DSD-855-A	Oral Fluid Multi-Drug 5 Drugs Test A Cassette	Oral Fluid	Cassette	25 T	/	CE	
DSD-865-A	Oral Fluid Multi-Drug 6 Drugs Test A Cassette	Oral Fluid	Cassette	25 T	/	CE	
DSD-875-A	Oral Fluid Multi-Drug 7 Drugs Test A Cassette	Oral Fluid	Cassette	25 T	/	CE	
DSD-885-A	Oral Fluid Multi-Drug 8 Drugs Test A Cassette	Oral Fluid	Cassette	25 T	/	CE	
DSD-895-A	Oral Fluid Multi-Drug 9 Drugs Test A Cassette	Oral Fluid	Cassette	25 T	/	CE	
DSD-827	Multi-Drug 2 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-837	Multi-Drug 3 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-847	Multi-Drug 4 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-857	Multi-Drug 5 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-867	Multi-Drug 6 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-877	Multi-Drug 7 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-887	Multi-Drug 8 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-897	Multi-Drug 9 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-8107	Multi-Drug 10 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-8117	Multi-Drug 11 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DSD-8127	Multi-Drug 12 Drugs Rapid Test Cup	Oral Fluid	Cup	20 T	/	CE	
DAC-X14	Acetaminophen (ACE) Rapid Test Panel	Powder	Panel	25 T	5000 ng/mL	CE	
DAM-X14	Amphetamine (AMP) Rapid Test Panel	Powder	Panel	25 T	1000 ng/mL	CE	
DAP-X14	$\alpha$ -Pyrrolidinovalerophenone ( $\alpha$ -PVP ) Rapid Test Panel	Powder	Panel	25 T	1000 ng/mL	CE	
DBA-X14	Barbiturate (BAR) Rapid Test Panel	Powder	Panel	25 T	300 ng/mL	CE	
DBU-X14	Buprenorphine (BUP) Rapid Test Panel	Powder	Panel	25 T	10 ng/mL	CE	
DBZ-X14	Benzodiazepines (BZO) Rapid Test Panel	Powder	Panel	25 T	300 ng/mL	CE	
DCFY-X14	Carfentanyl(CFYL) Rapid Test Panel	Powder	Panel	25 T	500 ng/mL	CE	
DCO-X14	COCAINE (COC) Rapid Test Panel	Powder	Panel	25 T	300 ng/mL	CE	
DCT-X14	Cotinine (COT) Rapid Test Panel	Powder	Panel	25 T	200 ng/mL	CE	
DDI-X14	Diazepam (DIA) Rapid Test Panel	Powder	Panel	25 T	300 ng/mL	CE	
DFY-X14	Fentanyl (FYL) Rapid Test Panel	Powder	Panel	25 T	200 ng/mL	CE	
DKE-X14	Ketamine (KET)Rapid Test Panel	Powder	Panel	25 T	1000 ng/mL	CE	
DMA-X14	6-Monoacetylmorphine (6-MAM)Rapid Test Panel	Powder	Panel	25 T	10 ng/mL	CE	
DMD-X14	Ecstasy (MDMA) Rapid Test Panel	Powder	Panel	25 T	500 ng/mL	CE	
DMDA-X14	Tenamfetamine (MDA) Rapid Test Panel	Powder	Panel	25 T	500 ng/mL	CE	
DMDP-X14	Methylenedioxypropyvalerone (MDPV) Rapid Test Panel	Powder	Panel	25 T	1000 ng/mL	CE	
DMP-X14	Methylphenidate(MPD) Rapid Test Panel	Powder	Panel	25 T	300 ng/mL	CE	

Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
DME-X14	Methamphetamine (MET) Rapid Test Panel	Powder	Panel	25 T	1000 ng/mL	CE	
DMO-X14	Morphine (MOP) Rapid Test Panel	Powder	Panel	25 T	300 ng/mL	CE	
DMQ-X14	Methaqualone (MQL) Rapid Test Panel	Powder	Panel	25 T	300 ng/mL	CE	
DMT-X14	Methadone (MTD) Rapid Test Panel	Powder	Panel	25 T	300 ng/mL	CE	
DOX-X14	Oxycodone (OXY) Rapid Test Panel	Powder	Panel	25 T	100 ng/mL	CE	
DOP-X14	Opiate (OPI) Rapid Test Panel	Powder	Panel	25 T	2000 ng/mL	CE	
DPC-X14	Phencyclidine (PCP) PCP Rapid Test Panel	Powder	Panel	25 T	25 ng/mL	CE	
DPP-X14	Propoxyphene (PPX) Rapid Test Panel	Powder	Panel	25 T	300 ng/mL	CE	
DSM-X14	Synthetic Marijuana (K2) Rapid Test Panel	Powder	Panel	25 T	50 ng/mL	CE	
DTC-X14	Tricyclic Antidepressants (TCA) Rapid Test Panel	Powder	Panel	25 T	1000 ng/mL	CE	
DTH-X14	Marijuana (THC) Rapid Test Panel	Powder	Panel	25 T	50 ng/mL	CE	
DTM-X14	Tramadol (TML) Rapid Test Panel	Powder	Panel	25 T	100 ng/mL	CE	
DLS-X14	Lysergic Acid Diethylamide (LSD) Rapid Test Panel	Powder	Panel	25 T	20 ng/mL	CE	
DZO-X14	Zolpidem(ZOL) Rapid Test Panel	Powder	Panel	25 T	50 ng/mL	CE	
DOA-X24	Multi-Drug 2 Drugs Rapid Test Panel	Powder	Panel	25 T	/	CE	
DOA-X34	Multi-Drug 3 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DOA-X44	Multi-Drug 4 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DOA-X54	Multi-Drug 5 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DOA-X64	Multi-Drug 6 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DOA-X74	Multi-Drug 7 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DOA-X84	Multi-Drug 8 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DOA-X94	Multi-Drug 9 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DOA-X104	Multi-Drug 10 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DOA-X114	Multi-Drug 11 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DOA-X124	Multi-Drug 12 Drugs Rapid Test Panel	Powder	Panel	25 T		CE	
DAM-402	AMP Rapid Test Cassette	WB/S/P	Cassette	25 T		80 ng/mL	CE
DBA-402	BAR Rapid Test Cassette	WB/S/P	Cassette	25 T	100 ng/mL	CE	
DBU-402	BUP Rapid Test Cassette	WB/S/P	Cassette	25 T	5 ng/mL	CE	
DBZ-402	BZO Rapid Test Cassette	WB/S/P	Cassette	25 T	100 ng/mL	CE	
DCO-402	COC Rapid Test Cassette	WB/S/P	Cassette	25 T	50 ng/mL	CE	
DCT-402	COT Rapid Test Cassette	WB/S/P	Cassette	25 T	10 ng/mL*,100 ng/mL	CE	
DED-402	EDDP Rapid Test Cassette	WB/S/P	Cassette	25 T	50 ng/mL	CE	
DMD-402	MDMA Rapid Test Cassette	WB/S/P	Cassette	25 T	50 ng/mL	CE	
DME-402	MET Rapid Test Cassette	WB/S/P	Cassette	25 T	70 ng/mL	CE	
DMO-402	MOP Rapid Test Cassette	WB/S/P	Cassette	25 T	40 ng/mL	CE	
DMT-402	MTD Rapid Test Cassette	WB/S/P	Cassette	25 T	40 ng/mL	CE	
DOX-402	OXY Rapid Test Cassette	WB/S/P	Cassette	25 T	20 ng/mL	CE	



*Product List-Your Reliable Rapid Test Partner*

9500 7th Street Unit M Rancho Cucamonga , CA 91730

Tel: 1-909-466-6892 Fax: 1-909-466-6892 E-mail: selina.wu@acrobiotech.com

[www.acrobiotech.com](http://www.acrobiotech.com)

Cat. No.	Product Description	Specimen	Format	Kit Size	Cut-Off	Status	Picture
DPP-402	PPX Rapid Test Cassette	WB/S/P	Cassette	25 T	100 ng/mL	CE	
DTC-402	TCA Rapid Test Cassette	WB/S/P	Cassette	25 T	300 ng/mL	CE	
DTH-402	THC Rapid Test Cassette	WB/S/P	Cassette	25 T	35 ng/mL*,50 ng/mL	CE	
DTM-402	TML Rapid Test Cassette	WB/S/P	Cassette	25 T	50 ng/mL	CE	
DOA-425	Multi-drug 2 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-435	Multi-drug 3 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-445	Multi-drug 4 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-455	Multi-drug 5 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-465	Multi-drug 6 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-475	Multi-drug 7 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-485	Multi-drug 8 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-495	Multi-drug 9 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-4105	Multi-drug 10 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-4115	Multi-drug 11 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	
DOA-4125	Multi-drug 12 drugs Test Cassette	WB/S/P	Cassette	25 T		CE	